






## Product data sheet REF 73.42.30

11.12.2025

General data			
article number	73.42.30	prolongation tube 30cm ang	
Unit	piece	Customs code (HS code)	90189084
GMDN	62467	GTIN	04046826134181
Risk category	1r (EU MDR)	UMDNS-Code	16-270
Classification	6 (EU MDR)	Base UDI-DI	4046826;B;R002;02;ah;3V
Net weight (kg)	0,020	Operating instructions	GATDR00X
Latex free	yes	Preparation instructions	
 			
Medicon media	450.00.37 (catalogue ENT )		
Complete picture		Detail images	
			
size	30cm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4301 (X5CrNi18-10)		
Screw	1.4021 (X20Cr13)		
Others	1.4305 (X8CrNiS18-9)		