



## Product data sheet REF 73.43.10

11.12.2025


General data			
article number	73.43.10	Killian laring fcps tip Ø 5mm str	
Unit	piece	Customs code (HS code)	90189084
GMDN	62467	GTIN	04046826134297
Risk category	1r (EU MDR)	UMDNS-Code	16-270
Classification	6 (EU MDR)	Base UDI-DI	4046826;B;R002;02;ah;3V
Net weight (kg)	0,006	Operating instructions	GATDR00X
Latex free	yes	Preparation instructions	



Medicon media 450.00.37 (catalogue ENT )

## Complete picture

## Detail images



size	Ø 5mm	Coating	
------	-------	---------	--

## Storage conditions

No special requirements

## Materials

	Material number (DIN designation)
Main part	1.4021 (X20Cr13)
Others	1.4310 (X10CrNi18-8)