






## Product data sheet REF 73.46.57

11.12.2025

| General data  |                                   |                          |                         |
|---|-----------------------------------|--------------------------|-------------------------|
| article number  | 73.46.57                          | Jackson laryng fcps 27cm |                         |
| Unit  | piece                             | Customs code (HS code)   | 90189084                |
| GMDN  | 62467                             | GTIN                     | 04046826134808          |
| Risk category   | 1r (EU MDR)                       | UMDNS-Code               | 16-270                  |
| Classification  | 6 (EU MDR)                        | Base UDI-DI              | 4046826;B;R002;03;ah;48 |
| Net weight (kg)   | 0,060                             | Operating instructions   | GATDR00X                |
| Latex free  | yes                               | Preparation instructions | MPK12-1                 |
|   |                                   |                          |                         |
| Medicon media   | 450.00.37 (catalogue ENT )        |                          |                         |
| Complete picture  |                                   | Detail images            |                         |
|   |                                   |                          |                         |
| size  | 27cm                              | Coating                  |                         |
| Storage conditions  |                                   |                          |                         |
| No special requirements   |                                   |                          |                         |
| Materials   |                                   |                          |                         |
|   | Material number (DIN designation) |                          |                         |
| Main part   | 1.4024 (X15Cr13)                  |                          |                         |