







Product data sheet REF 73.47.14

11.12.2025

General data			
article number	73.47.14	laryng mirror fig 2 Ø 14mm p/5	
Unit	box/5	Customs code (HS code)	90189084
GMDN	33431	GTIN	04046826134877
Risk category	1 (EU MDR)	UMDNS-Code	12-549
Classification	5 (EU MDR)	Base UDI-DI	4046826;N004;01;XD
Net weight (kg)	0,043	Operating instructions	GATDN00X
Latex free	yes	Preparation instructions	
 			
Medicon media	450.00.03 (general catalogue the art of surgery) 450.00.37 (catalogue ENT)		
Complete picture		Detail images	
 			
size	Ø 14mm	Coating	Vapor deposited coating on laryngeal mirrors
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4301 (X5CrNi18-10)		
Others	GLAS (Glas)		