







## Product data sheet REF 73.48.28

11.12.2025

General data			
article number	73.48.28	laryng mirror 22cm fig 9 Ø 28mm	
Unit	piece	Customs code (HS code)	90189084
GMDN	33431	GTIN	04046826135096
Risk category	1 (EU MDR)	UMDNS-Code	12-549
Classification	5 (EU MDR)	Base UDI-DI	4046826;N004;01;XD
Net weight (kg)	0,031	Operating instructions	GATDN00X
Latex free	yes	Preparation instructions	
 			
Medicon media	450.00.03 (general catalogue the art of surgery) 450.00.37 (catalogue ENT )		
Complete picture		Detail images	
 			
size	22cm Ø 28mm	Coating	Vapor deposited coating on laryngeal mirrors
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4301 (X5CrNi18-10)		
Others	GLAS (Glas)		