







Product data sheet REF 75.23.81

11.12.2025

General data			
article number	75.23.81	screwdriver 1.5mm 20cm	
Unit	piece	Customs code (HS code)	90189084
GMDN	33968	GTIN	04046826140120
Risk category	1r (EU MDR)	UMDNS-Code	13-517
Classification	6 (EU MDR)	Base UDI-DI	4046826;B;R005;03;ak;6A
Net weight (kg)	0,045	Operating instructions	GATDR00X
Latex free	yes	Preparation instructions	
 			
Medicon media	450.00.63 (catalogue Orthopaedics)		
Complete picture		Detail images	
 			
size	1.5mm 20cm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part Handle	1.4112 (X90CrMoV18) PPSU (Polyphenylsulfon)		