







Product data sheet REF 76.01.03

11.12.2025

General data			
article number	76.01.03	sut wire soft 0.30mm coil/10m	
Unit	piece	Customs code (HS code)	90211090
GMDN	56643	GTIN	04046826146542
Risk category	2b (EU MDD)	UMDNS-Code	16-104
Classification	8 (EU MDD)	Base UDI-DI	
Net weight (kg)	0,017	Operating instructions	G7080000
Latex free	yes	Preparation instructions	MPK01
  			
Medicon media	450.00.03 (general catalogue the art of surgery) 450.00.53 (catalogue cardio) 450.00.60 (catalogue plastic surgery) 450.00.62 (catalogue Medicon instruments Oral and Maxillofacial Surgery) 450.00.63 (catalogue Orthopaedics) 451.09.02 (brochure bimaxillary surgery)		
Complete picture		Detail images	
			
size	0.30mm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4441 (X2CrNiMo18-15-3)		