








Product data sheet REF 76.05.24

11.12.2025

| General data | | | |
|---|---|--------------------------|-------------------------|
| article number | 76.05.24 | wire guide 24cm | |
| Unit | piece | Customs code (HS code) | 90189084 |
| GMDN | 32864 | GTIN | 04046826146702 |
| Risk category | 1r (EU MDR) | UMDNS-Code | 11-925 |
| Classification | 6 (EU MDR) | Base UDI-DI | 4046826;B;R005;03;ak;6A |
| Net weight (kg) | 0,062 | Operating instructions | GATDR00X |
| Latex free | yes | Preparation instructions | |
|   | | | |
| Medicon media | 450.00.03 (general catalogue the art of surgery) 450.00.63 (catalogue Orthopaedics) | | |
| Complete picture | | Detail images | |
|    | | | |
| size | 24cm | Coating | |
| Storage conditions | | | |
| No special requirements | | | |
| Materials | | | |
| | Material number (DIN designation) | | |
| Main part | 1.4301 (X5CrNi18-10) | | |