



Product data sheet REF 76.06.14

25.08.2026

General data			
article number	76.06.14	Kirschner wire Ø 1.4mm 6cm tr/flat p/10	
Unit	box/10	Customs code (HS code)	90211090
GMDN	32854	GTIN	04046826146818
Risk category	2b (EU MDD)	UMDNS-Code	16-104
Classification	8 (EU MDD)	Base UDI-DI	4046826;B001;05;PL
Net weight (kg)	0,001	Operating instructions	G7080000
Latex free	yes	Preparation instructions	MPK01
			
Medicon media	450.00.03 (general catalogue the art of surgery) 450.00.63 (catalogue Orthopaedics)		
Complete picture		Detail images	
size	Ø 1.4mm 6cm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4441 (X2CrNiMo18-15-3)		