







## Product data sheet REF 76.07.80

11.12.2025

| General data  |   |  |                |
|---|---|--|----------------|
| article number  | 76.07.80  | Kirschner wire Ø 1.0mm 10cm tr/tr p/10 |                |
| Unit  | box/10  | Customs code (HS code)                 | 90211090       |
| GMDN  | 56643   | GTIN                                   | 04046826147525 |
| Risk category   | 2b (EU MDD)   | UMDNS-Code                             | 16-104         |
| Classification  | 8 (EU MDD)  | Base UDI-DI                            |                |
| Net weight (kg)   | 0,017   | Operating instructions                 | G7080000       |
| Latex free  | yes   | Preparation instructions               | MPK01          |
|    |   |  |                |
| Medicon media   | 450.00.03 (general catalogue the art of surgery)<br>450.00.63 (catalogue Orthopaedics ) |  |                |
| Complete picture  |   | Detail images                          |                |
|   |   |  |                |
| size  | Ø 1.0mm 10cm  | Coating                                |                |
| Storage conditions  |   |  |                |
| No special requirements   |   |  |                |
| Materials   |   |  |                |
|   | Material number (DIN designation)   |  |                |
| Main part   | 1.4441 (X2CrNiMo18-15-3)  |  |                |