







Product data sheet REF 76.11.58

11.12.2025

General data			
article number	76.11.58	Kirschner wire Ø 0.8mm 23cm tr/rd p/10	
Unit	box/10	Customs code (HS code)	90211090
GMDN	56643	GTIN	04046826149062
Risk category	2b (EU MDD)	UMDNS-Code	16-104
Classification	8 (EU MDD)	Base UDI-DI	
Net weight (kg)	0,025	Operating instructions	G7080000
Latex free	yes	Preparation instructions	MPK01
  			
Medicon media	450.00.03 (general catalogue the art of surgery) 450.00.63 (catalogue Orthopaedics)		
Complete picture		Detail images	
			
size	Ø 0.8mm 23cm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4441 (X2CrNiMo18-15-3)		