



Product data sheet REF 76.11.70

11.12.2025

General data			
article number	76.11.70	Kirschner wire Ø 2.0mm 23cm tr/rd p/10	
Unit	box/10	Customs code (HS code)	90211090
GMDN	56643	GTIN	04046826149154
Risk category	2b (EU MDD)	UMDNS-Code	16-104
Classification	8 (EU MDD)	Base UDI-DI	
Net weight (kg)	0,068	Operating instructions	G7080000
Latex free	yes	Preparation instructions	MPK01



Medicon media	450.00.03 (general catalogue the art of surgery) 450.00.63 (catalogue Orthopaedics)
---------------	---

Complete picture

Detail images

—

medicon®

size	Ø 2.0mm 23cm	Coating	
------	--------------	---------	--

Storage conditions

No special requirements

Materials

	Material number (DIN designation)
Main part	1.4441 (X2CrNiMo18-15-3)