







Product data sheet REF 76.15.18

11.12.2025

| General data | | | |
|---|---|--|----------------|
| article number | 76.15.18 | Kirschner wire Ø 1.8mm 31cm lancet/flat p/10 | |
| Unit | box/10 | Customs code (HS code) | 90211090 |
| GMDN | 56643 | GTIN | 04046826149673 |
| Risk category | 2b (EU MDD) | UMDNS-Code | 16-104 |
| Classification | 8 (EU MDD) | Base UDI-DI | |
| Net weight (kg) | 0,076 | Operating instructions | G7080000 |
| Latex free | yes | Preparation instructions | MPK01 |
|    | | | |
| Medicon media | 450.00.03 (general catalogue the art of surgery) 450.00.63 (catalogue Orthopaedics) | | |
| Complete picture | | Detail images | |
|  | | | |
| size | Ø 1.8mm 31cm | Coating | |
| Storage conditions | | | |
| No special requirements | | | |
| Materials | | | |
| | Material number (DIN designation) | | |
| Main part | 1.4441 (X2CrNiMo18-15-3) | | |