







## Product data sheet REF 76.15.22

11.12.2025

General data			
article number	76.15.22	Kirschner wire Ø 2.2mm 31cm lancet/flaz p/10	
Unit	box/10	Customs code (HS code)	90211090
GMDN	56643	GTIN	04046826149703
Risk category	2b (EU MDD)	UMDNS-Code	16-104
Classification	8 (EU MDD)	Base UDI-DI	
Net weight (kg)	0,105	Operating instructions	G7080000
Latex free	yes	Preparation instructions	MPK01
  			
Medicon media	450.00.03 (general catalogue the art of surgery) 450.00.63 (catalogue Orthopaedics )		
Complete picture		Detail images	
			
size	Ø 2.2mm 31cm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4441 (X2CrNiMo18-15-3)		