





Product data sheet REF 76.16.01

11.12.2025

General data			
article number	76.16.01	wire dispenser f 16cm Ø 1.8/2.0/2.2/2.5mm	
Unit	piece	Customs code (HS code)	90189084
GMDN	12143	GTIN	04046826149741
Risk category	1 (EU MDR)	UMDNS-Code	16-767
Classification	1 (EU MDR)	Base UDI-DI	4046826;N003;01;X2
Net weight (kg)	0,120	Operating instructions	GATDN00X
Latex free	yes	Preparation instructions	
 			
Medicon media			
Complete picture		Detail images	
size	16cm Ø 1.8/2.0/2.2/2.5mm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4301 (X5CrNi18-10)		
Screw	1.4305 (X8CrNiS18-9)		
Handle	1.4310 (X10CrNi18-8)		
Others	POM (Polyoxymethylen_(Polyacetal/_Polyformaldehyd)		