







## Product data sheet REF 76.16.03

11.12.2025

General data			
article number	76.16.03	Kirschner wire tube 16cm Ø 12mm	
Unit	piece	Customs code (HS code)	90189084
GMDN	12143	GTIN	04046826149758
Risk category	1 (EU MDR)	UMDNS-Code	16-767
Classification	1 (EU MDR)	Base UDI-DI	4046826;N003;03;X8
Net weight (kg)	0,100	Operating instructions	GATDN00X
Latex free	yes	Preparation instructions	
 			
Medicon media	450.00.03 (general catalogue the art of surgery) 450.00.63 (catalogue Orthopaedics )		
Complete picture		Detail images	
 			
size	16cm Ø 12mm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part Screw	1.4301 (X5CrNi18-10) 1.4305 (X8CrNiS18-9)		