







## Product data sheet REF 76.19.51

11.12.2025

General data			
article number	76.19.51	Steinmann nail tr/triang Ø 5.0mm 15cm center thread	
Unit	piece	Customs code (HS code)	90211090
GMDN	56643	GTIN	04046826150495
Risk category	2b (EU MDD)	UMDNS-Code	16-078
Classification	8 (EU MDD)	Base UDI-DI	
Net weight (kg)	0,020	Operating instructions	G7080000
Latex free	yes	Preparation instructions	MPK01
  			
Medicon media	450.00.03 (general catalogue the art of surgery) 450.00.63 (catalogue Orthopaedics )		
Complete picture		Detail images	
			
size	Ø 5.0mm 15cm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4441 (X2CrNiMo18-15-3)		