



## Product data sheet REF 76.21.90

11.12.2025

General data			
article number	76.21.90	Steinmann nail tr/tr Ø 4.0mm 23cm	
Unit	piece	Customs code (HS code)	90211090
GMDN	56643	GTIN	04046826150716
Risk category	2b (EU MDD)	UMDNS-Code	16-078
Classification	8 (EU MDD)	Base UDI-DI	
Net weight (kg)	0,022	Operating instructions	G7080000
Latex free	yes	Preparation instructions	MPK01



Medicon media	450.00.03 (general catalogue the art of surgery) 450.00.63 (catalogue Orthopaedics )
---------------	---

## Complete picture

## Detail images



size	Ø 4.0mm 23cm	Coating	
------	--------------	---------	--

## Storage conditions

No special requirements

## Materials

	Material number (DIN designation)
Main part	1.4441 (X2CrNiMo18-15-3)