



Product data sheet REF 76.23.95

11.12.2025

General data			
article number	76.23.95	Steinmann nail tr/triang Ø 4.5mm 23cm	
Unit	piece	Customs code (HS code)	90211090
GMDN	56643	GTIN	04046826151003
Risk category	2b (EU MDD)	UMDNS-Code	16-078
Classification	8 (EU MDD)	Base UDI-DI	
Net weight (kg)	0,020	Operating instructions	G7080000
Latex free	yes	Preparation instructions	MPK01
			
Medicon media	450.00.03 (general catalogue the art of surgery) 450.00.63 (catalogue Orthopaedics)		
Complete picture		Detail images	
			
size	Ø 4.5mm 23cm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4441 (X2CrNiMo18-15-3)		