



## Product data sheet REF 76.41.20

11.12.2025

| General data    |            |                                 |                    |
|-----------------|------------|---------------------------------|--------------------|
| article number  | 76.41.20   | tightening wrench f 76.41.01-07 |                    |
| Unit            | piece      | Customs code (HS code)          | 90189084           |
| GMDN            | 46525      | GTIN                            | 04046826151164     |
| Risk category   | 1 (EU MDR) | UMDNS-Code                      | 14-101             |
| Classification  | 1 (EU MDR) | Base UDI-DI                     | 4046826;N001;06;WT |
| Net weight (kg) | 0,275      | Operating instructions          | GATDN00X           |
| Latex free      | yes        | Preparation instructions        |                    |



Medicon media 450.00.03 (general catalogue the art of surgery)  
450.00.63 (catalogue Orthopaedics )

## Complete picture

## Detail images



## Storage conditions

No special requirements

## Materials

|           | Material number (DIN designation) |
|-----------|-----------------------------------|
| Main part | 1.4021 (X20Cr13)                  |