



Product data sheet REF 76.44.25

11.12.2025

General data			
article number	76.44.25	Steinmann pin 25cm	
Unit	piece	Customs code (HS code)	90211090
GMDN	56643	GTIN	04046826151348
Risk category	2b (EU MDD)	UMDNS-Code	16-078
Classification	8 (EU MDD)	Base UDI-DI	
Net weight (kg)	0,038	Operating instructions	G7080000
Latex free	yes	Preparation instructions	MPK01
  			
Medicon media			
Complete picture		Detail images	
size	25cm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4441 (X2CrNiMo18-15-3)		