






Product data sheet REF 76.45.73

11.12.2025

General data			
article number	76.45.73	wire cutting pliers 23cm	
Unit	piece	Customs code (HS code)	90189084
GMDN	32885	GTIN	04046826151522
Risk category	1 (EU MDR)	UMDNS-Code	13-041
Classification	1 (EU MDR)	Base UDI-DI	4046826;N001;03;WJ
Net weight (kg)	0,409	Operating instructions	GATDN00X
Latex free	yes	Preparation instructions	
 			
Medicon media	450.00.03 (general catalogue the art of surgery) 450.00.53 (catalogue cardio) 450.00.63 (catalogue Orthopaedics) 451.03.23 (brochure medicut scrs)		
Complete picture		Detail images	
			
size	23cm	Coating	Artificially created galvanic gold layer (max 0.04mm) on parts of the surface of instruments
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4111 (X110CrMoV15)		
Screw	1.4305 (X8CrNiS18-9)		
Others	1.4021 (X20Cr13)		