



Product data sheet REF 76.45.85

11.12.2025

General data			
article number	76.45.85	Easy Cut HM wire a pin cutter 20cm	
Unit	piece	Customs code (HS code)	90189084
GMDN	32885	GTIN	04046826151584
Risk category	1 (EU MDR)	UMDNS-Code	13-041
Classification	6 (EU MDR)	Base UDI-DI	4046826;N001;03;WJ
Net weight (kg)	0,350	Operating instructions	GATDN00X
Latex free	yes	Preparation instructions	
 			
Medicon media	450.00.03 (general catalogue the art of surgery) 450.00.53 (catalogue cardio) 450.00.63 (catalogue Orthopaedics) 451.03.23 (brochure medicut scrs)		
Complete picture		Detail images	
  			
size	20cm	Coating	Artificially created PVD layer of titanium nitrite e.g. Blackline scissors
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4021 (X20Cr13)		
Screw	1.4305 (X8CrNiS18-9)		
Others	1.4125 (X110CrMo17)		