






## Product data sheet REF 76.50.12

11.12.2025

General data			
article number	76.50.12	twist drill fig 2 Ø 3.0mm	
Unit	piece	Customs code (HS code)	90189084
GMDN	32390	GTIN	04046826151683
Risk category	1r (EU MDR)	UMDNS-Code	11-331
Classification	6 (EU MDR)	Base UDI-DI	4046826;B;R003;05;ai;5L
Net weight (kg)	0,007	Operating instructions	GATDR00X
Latex free	yes	Preparation instructions	MPK01-1
 			
Medicon media	450.00.03 (general catalogue the art of surgery) 450.00.63 (catalogue Orthopaedics )		
Complete picture		Detail images	
			
size	Ø 3.0mm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part Screw	1.4021 (X20Cr13) 1.4305 (X8CrNiS18-9)		