






## Product data sheet REF 76.69.20

11.12.2025

General data			
article number	76.69.20	twist drill Ø 2.0mm 8.5cm	
Unit	piece	Customs code (HS code)	90189084
GMDN	32390	GTIN	04046826152765
Risk category	1r (EU MDR)	UMDNS-Code	11-331
Classification	6 (EU MDR)	Base UDI-DI	4046826;B;R003;05;ai;5L
Net weight (kg)	0,003	Operating instructions	GATDR00X
Latex free	yes	Preparation instructions	MPK01-1
 			
Medicon media	450.00.03 (general catalogue the art of surgery) 450.00.63 (catalogue Orthopaedics )		
Complete picture		Detail images	
			
size	Ø 2.0mm 8.5cm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4021 (X20Cr13)		