



Product data sheet REF 76.69.45

11.12.2025

General data			
article number	76.69.45	twist drill Ø 4.5mm 13.0cm	
Unit	piece	Customs code (HS code)	90189084
GMDN	32390	GTIN	04046826152826
Risk category	1r (EU MDR)	UMDNS-Code	11-331
Classification	6 (EU MDR)	Base UDI-DI	4046826;B;R003;05;ai;5L
Net weight (kg)	0,010	Operating instructions	GATDR00X
Latex free	yes	Preparation instructions	MPK01-1



Medicon media	450.00.03 (general catalogue the art of surgery) 450.00.63 (catalogue Orthopaedics)
---------------	---

Complete picture

Detail images



size	Ø 4.5mm 13.0cm	Coating	
------	----------------	---------	--

Storage conditions

No special requirements

Materials

	Material number (DIN designation)
Main part	1.4034 (X46Cr13)