








Product data sheet REF 77.10.60

11.12.2025

General data			
article number	77.10.60	Lange Baby bone lever 17mm 17cm	
Unit	piece	Customs code (HS code)	90189084
GMDN	35048	GTIN	04046826153885
Risk category	1r (EU MDR)	UMDNS-Code	11-508
Classification	6 (EU MDR)	Base UDI-DI	4046826;B;R004;03;aj;5L
Net weight (kg)	0,045	Operating instructions	GATDR00X
Latex free	yes	Preparation instructions	MPK12-1
 			
Medicon media	450.00.03 (general catalogue the art of surgery) 450.00.63 (catalogue Orthopaedics)		
Complete picture		Detail images	
  			
size	17mm 17cm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4021 (X20Cr13)		