



## Product data sheet REF 77.10.81

11.12.2025

General data			
article number	77.10.81	Lange Hohmann modif bone lever 30mm 27cm	
Unit	piece	Customs code (HS code)	90189084
GMDN	35048	GTIN	04046826153922
Risk category	1r (EU MDR)	UMDNS-Code	11-508
Classification	6 (EU MDR)	Base UDI-DI	4046826;B;R004;03;aj;5L
Net weight (kg)	0,230	Operating instructions	GATDR00X
Latex free	yes	Preparation instructions	MPK24-1



Medicon media 450.00.03 (general catalogue the art of surgery)  
450.00.63 (catalogue Orthopaedics )

## Complete picture

## Detail images



size	30mm 27cm	Coating	
------	-----------	---------	--

## Storage conditions

No special requirements

## Materials

	Material number (DIN designation)
Main part	1.4021 (X20Cr13)