








Product data sheet REF 77.11.42

11.12.2025

General data			
article number	77.11.42	Verbrugge Hohmann bone lever 42mm 24cm	
Unit	piece	Customs code (HS code)	90189084
GMDN	35048	GTIN	04046826153991
Risk category	1r (EU MDR)	UMDNS-Code	11-508
Classification	6 (EU MDR)	Base UDI-DI	4046826;B;R004;03;aj;5L
Net weight (kg)	0,142	Operating instructions	GATDR00X
Latex free	yes	Preparation instructions	MPK20-1
 			
Medicon media	450.00.03 (general catalogue the art of surgery) 450.00.63 (catalogue Orthopaedics)		
Complete picture		Detail images	
  			
size	42mm 24cm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4021 (X20Cr13)		