



## Product data sheet REF 77.12.04

11.12.2025

General data			
article number	77.12.04	Bennett bone lever 45mm 24cm	
Unit	piece	Customs code (HS code)	90189084
GMDN	35048	GTIN	04046826154028
Risk category	1r (EU MDR)	UMDNS-Code	11-508
Classification	6 (EU MDR)	Base UDI-DI	4046826;B;R004;03;aj;5L
Net weight (kg)	0,175	Operating instructions	GATDR00X
Latex free	yes	Preparation instructions	MPK20-1



Medicon media 450.00.03 (general catalogue the art of surgery)  
450.00.63 (catalogue Orthopaedics )

## Complete picture

## Detail images



size	45mm 24cm	Coating	
------	-----------	---------	--

## Storage conditions

No special requirements

## Materials

	Material number (DIN designation)
Main part	1.4021 (X20Cr13)