




Product data sheet REF 77.60.31

11.12.2025

General data			
article number	77.60.31	mini endoscope Ø 1.9mm 60mm 30°	
Unit	piece	Customs code (HS code)	90189020
GMDN	34856	GTIN	04046826445157
Risk category	2a (EU MDD)	UMDNS-Code	10-198
Classification	6 (EU MDD)	Base UDI-DI	
Net weight (kg)	0,268	Operating instructions	G776000
Latex free	yes	Preparation instructions	
			
Medicon media			
Complete picture		Detail images	
size	Ø 1.9mm 60mm 30°	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4301 (X5CrNi18-10)		
Screw	1.4305 (X8CrNiS18-9)		
Handle	PEEK (Polyetheretherketon)		
Others	GLAS (Glas)		