






## Product data sheet REF 77.61.27

11.12.2025

General data			
article number	77.61.27	trocar triang Ø 4.0mm 175mm	
Unit	piece	Customs code (HS code)	90189084
GMDN	37144	GTIN	04046826154547
Risk category	1r (EU MDR)	UMDNS-Code	15-260
Classification	6 (EU MDR)	Base UDI-DI	4046826;B;R001;05;ag;48
Net weight (kg)	0,042	Operating instructions	GATDR00X
Latex free	yes	Preparation instructions	
 			
Medicon media	451.02.11 (brochure instruments f arthroscopy)		
Complete picture		Detail images	
			
size	Ø 4.0mm 175mm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4021 (X20Cr13)		