





Product data sheet REF 77.61.81

11.12.2025

| General data | | | |
|---|-----------------------------------|--|-------------------------|
| article number | 77.61.81 | obturator triang 35mm Ø 2mm f 77.61.80 | |
| Unit | piece | Customs code (HS code) | 90189084 |
| GMDN | 60909 | GTIN | 04046826450281 |
| Risk category | 1r (EU MDR) | UMDNS-Code | 15-260 |
| Classification | 6 (EU MDR) | Base UDI-DI | 4046826;B;R004;06;aj;6M |
| Net weight (kg) | 0,021 | Operating instructions | GATDR00X |
| Latex free | yes | Preparation instructions | |
|   | | | |
| Medicon media | | | |
| Complete picture | | Detail images | |
| | | | |
| size | 35mm Ø 2mm | Coating | |
| Storage conditions | | | |
| No special requirements | | | |
| Materials | | | |
| | Material number (DIN designation) | | |
| Main part | 1.4301 (X5CrNi18-10) | | |