







Product data sheet REF 77.89.19

11.12.2025

| General data | | | |
|---|--|------------------------------------|-------------------------|
| article number | 77.89.19 | punch Duck Bill Ø 3.4mm 45° lf cvd | |
| Unit | piece | Customs code (HS code) | 90189084 |
| GMDN | 35732 | GTIN | 04046826394165 |
| Risk category | 1r (EU MDR) | UMDNS-Code | 16-269 |
| Classification | 6 (EU MDR) | Base UDI-DI | 4046826;B;R003;06;ai;5X |
| Net weight (kg) | 0,063 | Operating instructions | GATDR00X |
| Latex free | yes | Preparation instructions | |
|   | | | |
| Medicon media | 451.02.11 (brochure instruments f arthroscopy) | | |
| Complete picture | | Detail images | |
|   | | | |
| size | Ø 3.4mm 45° | Coating | |
| Storage conditions | | | |
| No special requirements | | | |
| Materials | | | |
| | Material number (DIN designation) | | |
| Main part | 1.4034 (X46Cr13) | | |
| Screw | 1.4305 (X8CrNiS18-9) | | |
| Handle | 1.4021 (X20Cr13) | | |
| Others | 1.4301 (X5CrNi18-10) | | |