








## Product data sheet REF 77.95.19

11.12.2025

| General data   |  |                             |                         |
|--|--|-----------------------------|-------------------------|
| article number   | 77.95.19                                       | punch retrograde Ø 3.4mm lf |                         |
| Unit   | piece  | Customs code (HS code)      | 90189084                |
| GMDN   | 35732  | GTIN                        | 04046826394202          |
| Risk category  | 1r (EU MDR)                                    | UMDNS-Code                  | 16-269                  |
| Classification   | 6 (EU MDR)                                     | Base UDI-DI                 | 4046826;B;R003;06;ai;5X |
| Net weight (kg)  | 0,071  | Operating instructions      | GATDR00X                |
| Latex free   | yes  | Preparation instructions    |                         |
|    |  |                             |                         |
| Medicon media  | 451.02.11 (brochure instruments f arthroscopy) |                             |                         |
| Complete picture   |  | Detail images               |                         |
|    |  |                             |                         |
| size   | Ø 3.4mm  | Coating                     |                         |
| Storage conditions   |  |                             |                         |
| No special requirements  |  |                             |                         |
| Materials  |  |                             |                         |
|  | Material number (DIN designation)              |                             |                         |
| Main part  | 1.4034 (X46Cr13)                               |                             |                         |
| Screw  | 1.4305 (X8CrNiS18-9)                           |                             |                         |
| Handle   | 1.4021 (X20Cr13)                               |                             |                         |
| Others   | 1.4301 (X5CrNi18-10)                           |                             |                         |