



Product data sheet REF 78.03.01

11.12.2025

| General data | | | |
|-----------------|-------------|--------------------------|-------------------------|
| article number | 78.03.01 | bone reamer 16cm | |
| Unit | piece | Customs code (HS code) | 90189084 |
| GMDN | 45114 | GTIN | 04046826157357 |
| Risk category | 1r (EU MDR) | UMDNS-Code | 13-297 |
| Classification | 6 (EU MDR) | Base UDI-DI | 4046826;B;R003;05;ai;5L |
| Net weight (kg) | 0,095 | Operating instructions | GATDR00X |
| Latex free | yes | Preparation instructions | MPK12-1 |



| | |
|---------------|---|
| Medicon media | 450.00.03 (general catalogue the art of surgery) 450.00.63 (catalogue Orthopaedics) |
|---------------|---|

Complete picture

Detail images



| | | | |
|------|------|---------|--|
| size | 16cm | Coating | |
|------|------|---------|--|

Storage conditions

No special requirements

Materials

| | Material number (DIN designation) |
|-----------|-----------------------------------|
| Main part | 1.4021 (X20Cr13) |