








Product data sheet REF 78.12.13

11.12.2025

| General data | | | |
|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------|------------------------------------|-------------------------|
| article number | 78.12.13 | Steinmann tendon hold clamp 11.5cm | |
| Unit | piece | Customs code (HS code) | 90189084 |
| GMDN | 62467 | GTIN | 04046826157531 |
| Risk category | 1r (EU MDR) | UMDNS-Code | 15-221 |
| Classification | 6 (EU MDR) | Base UDI-DI | 4046826;B;R002;02;ah;3V |
| Net weight (kg) | 0,021 | Operating instructions | GATDR00X |
| Latex free | yes | Preparation instructions | |
|   | | | |
| Medicon media | 450.00.60 (catalogue plastic surgery) | | |
| Complete picture | | Detail images | |
|    | | | |
| size | 11.5cm | Coating | |
| Storage conditions | | | |
| No special requirements | | | |
| Materials | | | |
| | Material number (DIN designation) | | |
| Main part | 1.4021 (X20Cr13) | | |
| Screw | 1.4305 (X8CrNiS18-9) | | |
| Others | 1.4310 (X10CrNi18-8) | | |