








Product data sheet REF 78.12.15

11.12.2025

General data			
article number	78.12.15	tendon interlace clamp 15cm cvd	
Unit	piece	Customs code (HS code)	90189084
GMDN	62467	GTIN	04046826157548
Risk category	1r (EU MDR)	UMDNS-Code	15-221
Classification	6 (EU MDR)	Base UDI-DI	4046826;B;R002;02;ah;3V
Net weight (kg)	0,021	Operating instructions	GATDR00X
Latex free	yes	Preparation instructions	MPK23-1
 			
Medicon media	450.00.60 (catalogue plastic surgery)		
Complete picture		Detail images	
  			
size	15cm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4021 (X20Cr13)		