





## Product data sheet REF 78.15.04

11.12.2025

General data			
article number	78.15.04	Larsen tendon clamp 10cm 4mm	
Unit	piece	Customs code (HS code)	90189084
GMDN	62466	GTIN	04046826157630
Risk category	1r (EU MDR)	UMDNS-Code	15-221
Classification	6 (EU MDR)	Base UDI-DI	4046826;B;R002;03;ah;48
Net weight (kg)	0,011	Operating instructions	GATDR00X
Latex free	yes	Preparation instructions	
 			
Medicon media			
Complete picture		Detail images	
size	10cm 4mm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4024 (X15Cr13)		