






Product data sheet REF 78.22.50

11.12.2025

General data			
article number	78.22.50	Brand tendon stripper 50cm rd	
Unit	piece	Customs code (HS code)	90189084
GMDN	35380	GTIN	04046826157753
Risk category	1r (EU MDR)	UMDNS-Code	13-827
Classification	6 (EU MDR)	Base UDI-DI	4046826;B;R001;02;ag;37
Net weight (kg)	0,043	Operating instructions	GATDR00X
Latex free	yes	Preparation instructions	MPK12-1
 			
Medicon media	450.00.60 (catalogue plastic surgery)		
Complete picture		Detail images	
 			
size	50cm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4301 (X5CrNi18-10)		
Handle	2.0401 (CuZn39Pb3)		
Others	1.4310 (X10CrNi18-8)		