







Product data sheet REF 78.22.51

11.12.2025

General data			
article number	78.22.51	Brand tendon stripper 30cm rd	
Unit	piece	Customs code (HS code)	90189084
GMDN	35380	GTIN	04046826157760
Risk category	1r (EU MDR)	UMDNS-Code	13-827
Classification	6 (EU MDR)	Base UDI-DI	4046826;B;R001;02;ag;37
Net weight (kg)	0,024	Operating instructions	GATDR00X
Latex free	yes	Preparation instructions	MPK03-1
 			
Medicon media	450.00.60 (catalogue plastic surgery)		
Complete picture		Detail images	
 			
size	30cm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4301 (X5CrNi18-10)		
Handle	2.0401 (CuZn39Pb3)		
Others	1.4310 (X10CrNi18-8)		