








## Product data sheet REF 78.35.13

11.12.2025

| General data  |                                     |                          |                         |
|---|-------------------------------------|--------------------------|-------------------------|
| article number  | 78.35.13                            | bone lever 15mm 16cm     |                         |
| Unit  | piece                               | Customs code (HS code)   | 90189084                |
| GMDN  | 35048                               | GTIN                     | 04046826371609          |
| Risk category   | 1r (EU MDR)                         | UMDNS-Code               | 11-508                  |
| Classification  | 6 (EU MDR)                          | Base UDI-DI              | 4046826;B;R004;03;aj;5L |
| Net weight (kg)   | 0,000                               | Operating instructions   | GATDR00X                |
| Latex free  | yes                                 | Preparation instructions | MPK03-1                 |
|     |                                     |                          |                         |
| Medicon media   | 450.00.63 (catalogue Orthopaedics ) |                          |                         |
| Complete picture  |                                     | Detail images            |                         |
|    |                                     |                          |                         |
| size  | 15mm 16cm                           | Coating                  |                         |
| Storage conditions  |                                     |                          |                         |
| No special requirements   |                                     |                          |                         |
| Materials   |                                     |                          |                         |
|   | Material number (DIN designation)   |                          |                         |
| Main part   | 1.4021 (X20Cr13)                    |                          |                         |