



Product data sheet REF 79.68.88

11.12.2025

General data			
article number	79.68.88	oscillant electric autopsy saw 230V	
Unit	piece	Customs code (HS code)	90189084
GMDN	60909	GTIN	04046826311926
Risk category	KEIN MP (EU MDD)	UMDNS-Code	
Classification	NA (EU MDD)	Base UDI-DI	
Net weight (kg)	2,083	Operating instructions	G796888
Latex free	yes	Preparation instructions	KEIN MP
			
Medicon media	450.00.03 (general catalogue the art of surgery)		
Complete picture		Detail images	
			
size	230V	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	KEINE ANGABEN (Werkstoffangabe nicht vorhanden bzw. Stückliste)		