






## Product data sheet REF 87.30.92

11.12.2025

General data			
article number	87.30.92	monopolar cable 3m 4mm plug/plug protected	
Unit	piece	Customs code (HS code)	85446010
GMDN	60909	GTIN	04046826452773
Risk category	1 (EU MDR)	UMDNS-Code	11-493
Classification	1 (EU MDR)	Base UDI-DI	4046826;N002;06;X6
Net weight (kg)	0,000	Operating instructions	GATDN00X
Latex free	yes	Preparation instructions	
 			
Medicon media			
Complete picture		Detail images	
			
size	3m 4mm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	SILIKON (Poly(organo)siloxane)		