



## Product data sheet REF 87.30.95

11.12.2025

General data			
article number	87.30.95	monopolar cable 3m plug/socket	
Unit	piece	Customs code (HS code)	90189084
GMDN	35042	GTIN	04046826161453
Risk category	1 (EU MDR)	UMDNS-Code	11-493
Classification	1 (EU MDR)	Base UDI-DI	4046826;N002;06;X6
Net weight (kg)	0,084	Operating instructions	GATDN00X
Latex free	yes	Preparation instructions	



Medicon media	450.00.03 (general catalogue the art of surgery)
---------------	--

## Complete picture

## Detail images



size	3m	Coating	
------	----	---------	--

## Storage conditions

No special requirements

## Materials

	Material number (DIN designation)
Main part	SILIKON (Poly(organo)siloxane)
Handle	SILIKON (Poly(organo)siloxane)