






Product data sheet REF 87.33.51

11.12.2025

General data			
article number	87.33.51	electr extension hexagonal 17.5cm Ø 4.0mm	
Unit	piece	Customs code (HS code)	90189084
GMDN	44683	GTIN	04046826162139
Risk category	N.A. (EU MDD)	UMDNS-Code	
Classification	NA (EU MDD)	Base UDI-DI	
Net weight (kg)	0,005	Operating instructions	
Latex free	no	Preparation instructions	MPK02
			
Medicon media	450.00.03 (general catalogue the art of surgery) 450.00.37 (catalogue ENT)		
Complete picture		Detail images	
			
size	17.5cm Ø 4.0mm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4301 (X5CrNi18-10)		
Screw	2.0402 (CuZn40Pb2)		
Handle	PP_10%GF (Polypropylen_10%Glasfaser)		
Others	PEEK (Polyetheretherketon)		