








Product data sheet REF 88.00.20

11.12.2025

General data			
article number	88.00.20	bipolar cable Ø 12.5mm 3m	
Unit	piece	Customs code (HS code)	85446010
GMDN	35042	GTIN	04046826163167
Risk category	1 (EU MDR)	UMDNS-Code	11-493
Classification	1 (EU MDR)	Base UDI-DI	4046826;N002;06;X6
Net weight (kg)	0,090	Operating instructions	GATDN00X
Latex free	yes	Preparation instructions	
 			
Medicon media	450.00.03 (general catalogue the art of surgery) 450.00.37 (catalogue ENT) 450.00.40 (Katalog Nr 5 Dental) 450.00.46 (catalogue sinuscopy) 450.00.62 (catalogue Medicon instruments Oral and Maxillofacial Surgery) 451.02.51 (ENT basic sets Dr M W Winter Instr sets f ENT surgery) 451.06.28 (brochure bipolar non-sticking fcps)		
Complete picture		Detail images	
  			
size	Ø 12.5mm 3m	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	2.0060 (Cu-ETP)		
Handle	PA (Polyamid)		
Others	SI (Silicon_(Polysiloxane))		