








Product data sheet REF 88.00.21

11.12.2025

General data			
article number	88.00.21	bipolar cable Ø 12.5mm 5m	
Unit	piece	Customs code (HS code)	85446010
GMDN	35042	GTIN	04046826163174
Risk category	1 (EU MDR)	UMDNS-Code	11-493
Classification	1 (EU MDR)	Base UDI-DI	4046826;N002;06;X6
Net weight (kg)	0,127	Operating instructions	GATDN00X
Latex free	yes	Preparation instructions	
 			
Medicon media	450.00.03 (general catalogue the art of surgery) 450.00.37 (catalogue ENT) 450.00.62 (catalogue Medicon instruments Oral and Maxillofacial Surgery) 451.06.28 (brochure bipolar non-sticking fcps)		
Complete picture		Detail images	
  			
size	Ø 12.5mm 5m	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	2.0060 (Cu-ETP)		
Handle	PA (Polyamid)		
Others	SI (Silicon_(Polysiloxane))		