








## Product data sheet REF 88.00.23

11.12.2025

General data			
article number	88.00.23	bipolar cable Ø 8mm 3m	
Unit	piece	Customs code (HS code)	85446010
GMDN	35042	GTIN	04046826163198
Risk category	1 (EU MDR)	UMDNS-Code	11-493
Classification	1 (EU MDR)	Base UDI-DI	4046826;N002;06;X6
Net weight (kg)	0,086	Operating instructions	GATDN00X
Latex free	yes	Preparation instructions	
 			
Medicon media	450.00.03 (general catalogue the art of surgery) 450.00.37 (catalogue ENT ) 450.00.62 (catalogue Medicon instruments Oral and Maxillofacial Surgery) 451.06.28 (brochure bipolar non-sticking fcps)		
Complete picture		Detail images	
  			
size	Ø 8mm 3m	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	2.0060 (Cu-ETP)		
Others	PP_20%GF (Polypropylen_20%Glasfaser)		