



Product data sheet REF 88.00.29

11.12.2025

General data			
article number	88.00.29	cable bipolar Ø 4mm flat plug 5m 29mm	
Unit	piece	Customs code (HS code)	85446010
GMDN	35042	GTIN	04046826427443
Risk category	1 (EU MDR)	UMDNS-Code	11-493
Classification	1 (EU MDR)	Base UDI-DI	4046826;N002;06;X6
Net weight (kg)	0,000	Operating instructions	GATDN00X
Latex free	yes	Preparation instructions	



Medicon media	450.00.37 (catalogue ENT)
---------------	----------------------------

Complete picture

Detail images



size	Ø 4mm 5m	Coating	
------	----------	---------	--

Storage conditions

No special requirements

Materials

	Material number (DIN designation)
Main part	SI (Silicon_(Polysiloxane))
Others	PP_20%GF (Polypropylen_20%Glasfaser)